

Welcome to the era of
MOLECULAR ALLERGY for animals!



PAX

pet allergy xplorer



First quantitative macroarray
IgE test specifically designed
for animals

Over 200 allergen extracts
and molecular components


Better identification of allergen
cross-reactivities

Fully automated process, higher
level of standardisation

With CCD blocking and
2 blocking efficiency
detectors




Molecular Allergology: The future of IgE sensitisation detection



Molecular allergology is a state-of-the-art approach to the detection of sensitisations, whereby defined single allergen components are used for the determination of specific IgE in place of traditionally-used allergen extracts. The molecular components are recombinant proteins that provide a higher level of standardisation than allergen extracts and enable a more precise identification of IgE sensitisations. Molecular allergology tests are powerful tools that help pinpoint allergy triggers, thus facilitating risk assessment and therapy decisions.

Nextmune is bringing you the first molecular allergology platform for animals, the next-generation in allergen IgE serology:
PAX - Pet Allergy Xplorer




What are the main advantages of PAX?


- First quantitative multiplex macroarray specifically designed for companion animals
- Over 200 allergen extracts and components included = lower testing cost per allergen
- Fully automated process = higher level of standardization (same result if tested multiple times)
- With CCD blocking and 2 blocking efficiency detectors
- Only 0.5 ml of serum needed per test
- Expected increase in serological test sensitivity due to a higher concentration of molecular allergens
- Identification of "primary" sensitizing allergens
- Identification of allergen cross-reactivities
- Selection of relevant allergens for specific immunotherapy



Allergen
extract




Der f 2
NPC2 family




Der f 3
trypsin



Der f 10
tropomyosin



Der f 15
chitinase



Der f 1
cysteine protease

The PAX results are clearly set out, easy to interpret and include the following information:

- Summary of detectable sensitizations
- Interpretation summary and treatment recommendation
- Detailed results per extract and components
- Detailed interpretation with Information about allergenicity and relevance, time of the year, possible cross-reactivities and treatment indication for each allergen

PAX Complete result

The screenshots show the following sections of the PAX Complete result:

- Environmental:** A table listing various allergens such as Grass Pollen, Weeds, Dander & Epithelia, and Insects & Venoms, with their respective concentrations and allergenicity levels.
- Summary and immunotherapy recommendation:** A section providing a high-level overview of the patient's sensitization profile and recommending specific immunotherapy options (Option 1, Option 2, or a combination).
- Detailed allergen results:** A comprehensive list of allergens with detailed information on their allergenicity, relevance, and potential cross-reactivities.
- Interpretation - Support:** A section providing detailed clinical guidance, including treatment recommendations and information on how to interpret the results in the context of the patient's overall health and environment.

PAX Screening result

The screenshot shows the 'Screening Environmental' section of the PAX results. It displays a 'Positive' result for allergens, indicating that the patient has an elevated level of high-molecular-weight allergens. A 'PAX EXPAND' button is visible, suggesting the next step in the diagnostic process.



- Guaranteed 100% reliable screening test
- Fast results
- Continuous support and advice with our vet allergy experts

PAX Allergens: Components & Extracts

	Common name	Scientific name	Extracts & Components
Grass Pollens	Bermuda grass	<i>Cynodon dactylon</i>	Cyn d * rCyn d 1
	Orchard grass	<i>Dactylis glomerata</i>	Dac g *
	Meadow fescue	<i>Festuca pratensis</i>	Fes p *
	Perennial ryegrass	<i>Lolium perenne</i>	rLol p 1
	Timothy	<i>Phleum pratense</i>	rPhl p 1
			rPhl p 2
			rPhl p 5.0101
			rPhl p 6
			rPhl p 7
	Kentucky blue grass	<i>Poa pratensis</i>	Poa p *
Ryegrass, cultivated	<i>Secale cereale</i>	Sec c_pollen *	
Tree Pollens	Alder	<i>Alnus glutinosa</i>	Aln g *
			rAln g 1
			rAln g 4
	Silver birch	<i>Betula verrucosa</i>	Bet v *
			rBet v 1
			rBet v 2
			rBet v 6
	Hazel	<i>Corylus avellana</i>	Cor a_pollen * rCor a 1.0103
	Cypress	<i>Cupressus sempervirens</i>	Cup s *
	Beech	<i>Fagus sylvatica</i>	rFag s 1
	Ash	<i>Fraxinus excelsior</i>	Fra e * rFra e 1
	Privet	<i>Ligustrum vulgare</i>	Lig v *
	Olive tree	<i>Olea Europaea</i>	Ole e *
			nOle e 1
rOle e 7			
London plane tree	<i>Platanus acerifolia</i>	rPla a 1	
		nPla a 2 rPla a 3	
Cottonwood	<i>Populus nigra</i>	Pop n *	
Elm	<i>Ulmus campestris</i>	Ulm c *	
Weed Pollens	Ragweed	<i>Ambrosia artemisiifolia</i>	Amb a *
			rAmb a 1
			rAmb a 4
	Mugwort	<i>Artemisia vulgaris</i>	Art v *
			rArt v 1.0101
			rArt v 3.0201
	Lamb's quarter	<i>Chenopodium album</i>	Che a * rChe a 1
	Wall pellitory	<i>Parietaria judaica</i>	Par j * rPar j 2
	Ribwort / Plantain	<i>Plantago lanceolata</i>	Pla l * rPla l 1
	Dock/Sorrel	<i>Rumex crispus / acetosella</i>	Rum c / * Rum a
	Russian thistle	<i>Salsola kali</i>	Sal k * rSal k 1
	Nettle	<i>Urtica dioica</i>	Urt d *

	Common name	Scientific name	Extracts & Components
Danders & Epithelia	Cattle	<i>Bos domesticus</i>	rBos d 2
	Dog	<i>Canis familiaris</i>	rCan f 1
			rCan f 2
			nCan f 3
			rCan f 4
			rCan f 6
	Guinea pig	<i>Cavia porcellus</i>	rCav p 1
	Horse	<i>Equus caballus</i>	rEqu c 1
			nEqu c 3
			rEqu c 4
Cat	<i>Felis domesticu</i>	rFel d 1	
		nFel d 2	
		rFel d 4	
		rFel d 7	
Mouse	<i>Mus musculus</i>	rMus m 1	
Rabbit	<i>Oryctolagus cuniculus</i>	rOry c 1 rOry c 2 rOry c 3	
Mites & Cockroaches	Acarus siro	<i>Acarus siro</i>	Aca s *
	German cockroach	<i>Blattella germanica</i>	rBla g 1
			rBla g 2
			rBla g 4
			rBla g 5 rBla g 9
	Flea	<i>Ctenocephalides felis</i>	Cte f 1
	Dermatophagoides farinae	<i>Dermatophagoides farinae</i>	Der f *
			rDer f 1
			rDer f 2
			rDer f 15
			rDer f 18
			rDer f 18
	Dermatophagoides pteronyssinus	<i>Dermatophagoides pteronyssinus</i>	Der p *
			rDer p 1
			rDer p 2
			rDer p 5
			rDer p 7
			rDer p 10
rDer p 11			
rDer p 20			
rDer p 21			
rDer p 23			
Glycyphagus domesticus	<i>Glycyphagus domesticus</i>	rGly d 2	
Lepidoglyphus destructor	<i>Lepidoglyphus destructor</i>	Lep d *	
		rLep d 2	
Tyrophagus putrescentiae	<i>Tyrophagus putrescentiae</i>	Tyr p *	
		rTyr p 2	
Alternaria alternata	<i>Alternaria alternata</i>	Alt a *	
		rAlt a 1 rAlt a 6	

	Common name	Scientific name	Extracts & Components
Moulds & Yeasts	Aspergillus fumigatus	<i>Aspergillus fumigatus</i>	Asp f *
			rAsp f 1
			rAsp f 3
			rAsp f 4 rAsp f 6
	Cladosporium herbarum	<i>Cladosporium herbarum</i>	Cla h * rCla h 8
Malassezia pachydermatis	<i>Malassezia pachydermatis</i>	Mala p *	
Malassezia sympodialis	<i>Malassezia sympodialis</i>	rMala s 1	
		rMala s 9	
		rMala s 5	
		rMala s 6	
		rMala s 11	
Insect Venoms	Honey bee venom	<i>Apis mellifera</i>	Api m *
			nApi m 1
			Api m 2
			Api m 3 Api m 5 rApi m 10
	Long-headed wasp venom	<i>Dolichovespula spp.</i>	Dol spp *
Paper wasp venom	<i>Polistes dominulus</i>	Pol d * rPol d 5	
Fire ant venom	<i>Solenopsis richteri & Solenopsis invicta</i>	Sol spp *	
Common wasp venom	<i>Vespula vulgaris</i>	Ves v *	
		rVes v 1 rVes v 5	
Foods	Oat	<i>Avena sativa</i>	Ave s *
	Buckwheat	<i>Fagopyrum esculentum</i>	Fag e *
			nFag e 2
	Sunflower seed	<i>Helianthus annuus</i>	Hel a *
	Barley	<i>Hordeum vulgare</i>	Hor v *
	Rice	<i>Oryza sativa</i>	Ory s Ory s_GLUB1
	Millet	<i>Panicum miliaceum</i>	Pan m *
	Rye, cultivated	<i>Secale cereale</i>	Sec c_flour * Tri a *
	Wheat	<i>Triticum aestivum</i>	rTri a 14
			rTri a 19
			nTri a aA_TI
	Corn, cereal	<i>Zea mays</i>	Zea m *
rZea m 14 Zea m_GBSSI			
Apple	<i>Malus domestica</i>	rMal d 1	
		nMal d 2	
		rMal d 3	
Peanut	<i>Arachis hypogaea</i>	nAra h 1	
		rAra h 2	
		nAra h 3	
		rAra h 5	
		rAra h 6	
		rAra h 8 rAra h 9 rAra h 15	

* Extract

Common name	Scientific name	Extracts & Components
Soy	<i>Glycine max</i>	Gly m *
		rGly m 4
		rGly m 5
		nGly m 6
		nGly m 8
Lentil	<i>Lens culinaris</i>	Len c *
		Len c 1
		Len c 2
		Len c 3
Pea	<i>Pisum sativum</i>	Pis s *
		Pis s 1
		Pis s 2
		Pis s 3
Cow's milk	<i>Bos domesticus</i>	Bos d_milk *
		nBos d 4
		nBos d 5
		nBos d 8
Egg white	<i>Gallus domesticus</i>	Gal d_white *
		nGal d 1
		nGal d 2
		nGal d 3
Egg yolk	<i>Gallus domesticus</i>	Gal d_yolk *
		nGal d 5

Common name	Scientific name	Extracts & Components
Beef	<i>Bos domesticus</i>	Bos d_meat *
		nBos d 6
		Bos d 7
		Bos d_ACTA1
		Bos d_LDHA
Horse	<i>Equus caballus</i>	Equ c_meat *
Rabbit	<i>Oryctolagus spp.</i>	Ory_c_meat *
		Ory_c_CKM
		Ory_c_GAPDH
		Ory_c_PGM1
		Ory_c_PKM
Lamb	<i>Ovis aries</i>	Ovi a_meat *
		Ovi a_IgG
Pig	<i>Sus domesticus</i>	Sus d_meat *
Chicken	<i>Gallus domesticus</i>	rSus d 1
		Gal d_meat *
		Gal d 7
		Gal d 9
Turkey	<i>Meleagris gallopavo</i>	Gal d_PKM
		Mel g *
Mealworm	<i>Tenebrio molitor</i>	Ten m *
Herring, Atlantic	<i>Clupea harengus</i>	Clu h *
		rClu h 1

Common name	Scientific name	Extracts & Components
Cod, Atlantic	<i>Gadus morhua</i>	Gad m *
		nGad m 1
		Gad m 2+3
		Gad m 4
		Sal s *
Salmon, Atlantic	<i>Salmo salar</i>	Sal s 1
		Sal s 2
		Sal s 3
		Sal s 4
		Sal s 6
		Sal s 7
		Sal s 8
		Sal s 9
Mackerel, Atlantic	<i>Scomber scombrus</i>	Sco s *
		rSco s 1
Tuna	<i>Thunnus albacares</i>	Thu a *
Carrot	<i>Daucus carota</i>	Thu a 1
		Dau c *
Tomato	<i>Solanum lycopersicum</i>	rDau c 1
		Sola l *
Potato	<i>Solanum tuberosum</i>	rSola l 6
		Sol t *
		Sol t 2
		Sol t_GBSSI

* Extract



NextView is a newly developed portal where you can manage all your allergy samples, PAX results, immunotherapy orders, reorders, and much more.

With Nextview you can:

- Follow the status of your samples
- Access all your samples information, if they are in transit, being tested, and when results are expected.
- Easily find all results in one location
- Easily expand your screening results to complete panels
- Forward results directly to your customers for convenience
- Select and order a recommended treatment option with only one click
- Easily access your order history and re-order treatments with only one click
- Request a samples pick up (coming soon)
- Access your treatment reminder system (coming soon)
- Easily access your invoices (coming soon)



 **nextmune**



Nextmune | info.eu@nextmune.com | www.nextmune.com