

Welcome to the era of
MOLECULAR ALLERGY for animals!



PAX
pet allergy xplorer



First quantitative macroarray
IgE test specifically designed
for animals

Over 200 allergen extracts
and molecular components


Better identification of allergen
cross-reactivities

Fully automated process, higher
level of standardisation

With CCD blocking and
2 blocking efficiency
detectors




Molecular Allergology: The future of IgE sensitisation detection





Molecular allergology is a state-of-the-art approach to the detection of sensitisations, whereby defined single allergen components are used for the determination of specific IgE in place of traditionally-used allergen extracts. The molecular components are recombinant proteins that provide a higher level of standardisation than allergen extracts and enable a more precise identification of IgE sensitisations. Molecular allergology tests are powerful tools that help pinpoint allergy triggers, thus facilitating risk assessment and therapy decisions.

Nextmune is bringing you the first molecular allergology platform for animals, the next-generation in allergen-specific IgE serology.



What are the main advantages of PAX?

- First quantitative multiplex macroarray specifically designed for companion animals
- Over 200 allergen extracts and components included = lower testing cost per allergen
- Fully automated process = higher level of standardization (same result if tested multiple times)
- With CCD blocking and 2 blocking efficiency detectors
- Only 0.5 ml of serum needed
- Expected increase in serological test sensitivity due to a higher concentration of molecular allergens
- Identification of "primary" sensitizing allergens
- Identification of allergen cross-reactivities
- Selection of relevant allergens for specific immunotherapy



Molecular Allergology: The era of individual components

	Common name	Scientific name	Extracts & Components
Grass Pollens	Bermuda grass	<i>Cynodon dactylon</i>	Cyn d * rCyn d 1
	Orchard grass	<i>Dactylis glomerata</i>	Dac g *
	Meadow fescue	<i>Festuca pratensis</i>	Fes p *
	Perennial ryegrass	<i>Lolium perenne</i>	rLol p 1
	Timothy	<i>Phleum pratense</i>	rPhl p 1
			rPhl p 2
			rPhl p 5.0101
			rPhl p 6
Kentucky blue grass	<i>Poa pratensis</i>	rPhl p 7	
		rPhl p 12	
Ryegrass, cultivated	<i>Secale cereale</i>	Sec c_pollen *	
Tree Pollens	Alder	<i>Alnus glutinosa</i>	Aln g *
			rAln g 1
			rAln g 4
	Silver birch	<i>Betula verrucosa</i>	Bet v *
			rBet v 1
			rBet v 2
	Hazel	<i>Corylus avellana</i>	rBet v 6
			Cor a_pollen *
			rCor a 1.0103
	Cypress	<i>Cupressus sempervirens</i>	Cup s *
	Beech	<i>Fagus sylvatica</i>	rFag s 1
	Ash	<i>Fraxinus excelsior</i>	Fra e *
			rFra e 1
	Privet	<i>Ligustrum vulgare</i>	Lig v *
	Olive tree	<i>Olea Europaea</i>	Ole e *
nOle e 1			
rOle e 7			
London plane tree	<i>Platanus acerifolia</i>	rOle e 9	
		rPla a 1	
		nPla a 2	
Cottonwood	<i>Populus nigra</i>	rPla a 3	
		Pop n *	
Elm	<i>Ulmus campestris</i>	Ulm c *	
Weed Pollens	Ragweed	<i>Ambrosia artemisiifolia</i>	Amb a *
			rAmb a 1
			rAmb a 4
	Mugwort	<i>Artemisia vulgaris</i>	Art v *
			rArt v 1.0101
			rArt v 3.0201
	Lamb's quarter	<i>Chenopodium album</i>	Che a *
	Wall pellitory	<i>Parietaria judaica</i>	rChe a 1
			Par j *
	Ribwort / Plantain	<i>Plantago lanceolata</i>	rPar j 2
			Pla l *
	Dock/Sorrel	<i>Rumex crispus / acetosella</i>	rPla l 1
Rum c / *			
Russian thistle	<i>Salsola kali</i>	Rum a	
		Sal k *	
Nettle	<i>Urtica dioica</i>	rSal k 1	
		Urt d *	

	Common name	Scientific name	Extracts & Components
Danders & Epithelia	Cattle	<i>Bos domesticus</i>	rBos d 2
	Dog	<i>Canis familiaris</i>	rCan f 1
			rCan f 2
			nCan f 3
			rCan f 4
			rCan f 6
			Can f_maleurine (including Can f 5) *
	Guinea pig	<i>Cavia porcellus</i>	rCan f Fel d 1 like
			rCav p 1
	Horse	<i>Equus caballus</i>	rEqu c 1
nEqu c 3			
rEqu c 4			
Cat	<i>Felis domesticu</i>	rFel d 1	
		nFel d 2	
		rFel d 4	
		rFel d 7	
Mouse	<i>Mus musculus</i>	rMus m 1	
Rabbit	<i>Oryctolagus cuniculus</i>	rOry c 1	
		rOry c 2	
		rOry c 3	
Mites & Cockroaches	Acarus siro	<i>Acarus siro</i>	Aca s *
	German cockroach	<i>Blattella germanica</i>	rBla g 1
			rBla g 2
			rBla g 4
	Flea	<i>Ctenocephalides felis</i>	rBla g 5
			rBla g 9
			Cte f 1
	Dermatophagoides farinae	<i>Dermatophagoides farinae</i>	Der f *
			rDer f 1
			rDer f 2
	Dermatophagoides pteronyssinus	<i>Dermatophagoides pteronyssinus</i>	rDer f 15
			rDer f 18
			Der p *
			rDer p 1
			rDer p 2
			rDer p 5
			rDer p 7
			rDer p 10
rDer p 11			
rDer p 20			
rDer p 21			
rDer p 23			
Glycyphagus domesticus	<i>Glycyphagus domesticus</i>	rGly d 2	
Lepidoglyphus destructor	<i>Lepidoglyphus destructor</i>	Lep d *	
		rLep d 2	
Tyrophagus putrescentiae	<i>Tyrophagus putrescentiae</i>	Tyr p *	
		rTyr p 2	
Alternaria alternata	<i>Alternaria alternata</i>	Alt a *	
		rAlt a 1	
		rAlt a 6	

	Common name	Scientific name	Extracts & Components	
Moulds & Yeasts	Aspergillus fumigatus	<i>Aspergillus fumigatus</i>	Asp f *	
			rAsp f 1	
			rAsp f 3	
			rAsp f 4	
			rAsp f 6	
	Cladosporium herbarum	<i>Cladosporium herbarum</i>	Cla h *	
Malassezia pachydermatis	<i>Malassezia pachydermatis</i>	rCla h 8		
		Mala p *		
		rMala s 1		
		rMala s 9		
		rMala s 5		
Malassezia sympodialis	<i>Malassezia sympodialis</i>	rMala s 6		
		rMala s 11		
		Honey bee venom	<i>Apis mellifera</i>	Api m *
				nApi m 1
				Api m 2
				Api m 3
Api m 5				
rApi m 10				
Insect Venoms	Long-headed wasp venom	<i>Dolichovespula spp.</i>	Dol spp *	
	Paper wasp venom	<i>Polistes dominulus</i>	Pol d *	
	Fire ant venom	<i>Solenopsis richteri & Solenopsis invicta</i>	rPol d 5	
			Sol spp *	
	Common wasp venom	<i>Vespula vulgaris</i>	Ves v *	
			rVes v 1	
rVes v 5				
Foods	Oat	<i>Avena sativa</i>	Ave s *	
	Buckwheat	<i>Fagopyrum esculentum</i>	Fag e *	
			nFag e 2	
	Sunflower seed	<i>Helianthus annuus</i>	Hel a *	
	Barley	<i>Hordeum vulgare</i>	Hor v *	
	Rice	<i>Oryza sativa</i>	Ory s	
			Ory s_GLUB1	
	Millet	<i>Panicum miliaceum</i>	Pan m *	
	Rye, cultivated	<i>Secale cereale</i>	Sec c_flour *	
	Wheat	<i>Triticum aestivum</i>	Tri a *	
			rTri a 14	
			rTri a 19	
	Corn, cereal	<i>Zea mays</i>	nTri a aA_TI	
			Zea m *	
	Apple	<i>Malus domestica</i>	rZea m 14	
Zea m_GBSSI				
rMal d 1				
Peanut	<i>Arachis hypogaea</i>	rMal d 2		
		rMal d 3		
		nAra h 1		
Peanut	<i>Arachis hypogaea</i>	rAra h 2		
		rAra h 3		
		rAra h 5		
		rAra h 6		
		rAra h 8		
		rAra h 9		
rAra h 15				

* Extract

Common name	Scientific name	Extracts & Components
Soy	<i>Glycine max</i>	Gly m *
		rGly m 4
		rGly m 5
		nGly m 6
		nGly m 8
Lentil	<i>Lens culinaris</i>	Len c *
		Len c 1
		Len c 2
		Len c 3
Pea	<i>Pisum sativum</i>	Pis s *
		Pis s 1
		Pis s 2
		Pis s 3
Cow's milk	<i>Bos domesticus</i>	Bos d_milk *
		nBos d 4
		nBos d 5
		nBos d 8
Egg white	<i>Gallus domesticus</i>	Gal d_white *
		nGal d 1
		nGal d 2
		nGal d 3
		nGal d 4
Egg yolk	<i>Gallus domesticus</i>	Gal d_yolk * nGal d 5

Common name	Scientific name	Extracts & Components
Beef	<i>Bos domesticus</i>	Bos d_meat *
		nBos d 6
		Bos d 7
		Bos d_ACTA1
		Bos d_LDHA
Horse	<i>Equus caballus</i>	Equ c_meat *
Rabbit	<i>Oryctolagus spp.</i>	Ory c_meat *
		Ory c_CKM
		Ory c_GAPDH
		Ory c_PGM1
		Ory c_PKM
Lamb	<i>Ovis aries</i>	Ovi a_meat * Ovi a_lgG
Pig	<i>Sus domesticus</i>	Sus d_meat * rSus d 1
Chicken	<i>Gallus domesticus</i>	Gal d_meat *
		Gal d 7
		Gal d 9 Gal d_PKM
Turkey	<i>Meleagris gallopavo</i>	Mel g *
Mealworm	<i>Tenebrio molitor</i>	Ten m *
Herring, Atlantic	<i>Clupea harengus</i>	Clu h * rClu h 1

Common name	Scientific name	Extracts & Components
Cod, Atlantic	<i>Gadus morhua</i>	Gad m *
		nGad m 1
		Gad m 2+3
		Gad m 4
Salmon, Atlantic	<i>Salmo salar</i>	Sal s *
		Sal s 1
		Sal s 2
		Sal s 3
		Sal s 4
		Sal s 6
		Sal s 7
		Sal s 8
Mackerel, Atlantic	<i>Scomber scombrus</i>	Sco s *
		rSco s 1
Tuna	<i>Thunnus albacares</i>	Thu a * Thu a 1
Carrot	<i>Daucus carota</i>	Dau c * rDau c 1
Tomato	<i>Solanum lycopersicum</i>	Sola l * rSola l 6
Potato	<i>Solanum tuberosum</i>	Sol t *
		Sol t 2 Sol t_GBSSI

* Extract



NextView is a newly developed portal where you can manage all your allergy samples, PAX results, immunotherapy orders and reorders, and much more.

- With Nextview you can:
- Follow the status of your samples
- Access all your samples information, if they are in transit, being tested, and when results are expected.
- Easily find all results in one location
- Easily expand your screening results to complete panels
- Forward results directly to your customers for convenience
- Select and order a recommended treatment option with only one click
- Easily access your order history and re-order treatments with only one click
- Request a samples pick up (coming soon)
- Access your treatment reminder system (coming soon)
- Easily access your invoices (coming soon)



Register at www.nextmunelaboratories.co.uk

The PAX results are clearly set out, easy to interpret and include the following information:

- Summary of detectable sensitizations
- Interpretation summary and treatment recommendation
- Detailed results per extract and components
- Detailed interpretation with Information about allergenicity and relevance, time of the year, possible cross-reactivities and treatment indication for each allergen

In addition to all of this, we guarantee a 100% reliable screening test, fast results and continuous support and advice with our vet allergy experts.

PAX Complete result

The PAX Complete result report is divided into four main sections:

- Environmental Summary:** Provides a quick overview of detectable sensitizations, categorized by environmental allergens (Grass Pollen, Tree Pollen, Mites & Cockroaches), animal dander & epithelia, molds & yeast, and insects & venoms.
- Summary and Immunotherapy recommendation:** Offers a detailed interpretation of the results, including allergenicity and relevance information, and provides specific immunotherapy recommendations.
- Environmental Panel:** A comprehensive table listing all allergen extracts and their components, along with their respective IgE levels and clinical significance.
- Interpretation - Support:** Provides detailed clinical guidance, including advice on treatment options and further diagnostic steps.

PAX Screening result

The PAX Screening result report includes a section for environmental screening with the following details:

- Screening Environmental:** A section for allergen (extracts and components) screening.
- Positive for allergens when?** A question regarding the timing of allergen detection.
- Did we do an extra search to speed the result?** A question about additional testing procedures.
- PAX EXPAND:** An option to expand the screening results for more detailed information.

PAX - Pet Allergy Xplorer: The first molecular serum IgE specific test for pets

Dr. Thierry Olivry says:

Molecular allergology is a state-of-the-art approach to the detection of sensitisations, whereby defined single allergen components are used for the determination of specific IgE in place of traditionally-used allergen extracts. The molecular components are purified or recombinant proteins that provide a higher level of standardisation than allergen extracts and enable a more precise identification of IgE sensitisations. Pet Allergy Xplorer (PAX) is the first commercial serological IgE-specific test that uses allergen extracts and molecular components to identify which allergens are affecting pets.

THE PROBLEM

Traditionally, allergy testing in veterinary medicine involves placing an allergen extract on an enzyme linked immunoassay (ELISA) plate to incubate the serum and then administering a reagent that recognizes immunoglobulin E (IgE). The resultant color reaction indicates how much IgE is present. This technique, which is used by veterinary laboratories world-wide, has not changed for decades. However, results can vary considerably, depending on the extract used, and false negatives can occur if clinically relevant protein allergen concentrations are not sufficient. For example, the extract for the house dust mite, a common pet allergen, is made by grinding the mite, adding solvents to release the allergenic proteins, and purifying the proteins.

The house dust mite contains more than 10,000 proteins, but only about 40 that cause an allergic reaction are recognized. This means a low percentage of the allergy-causing proteins are seen when an extract is evaluated, especially if a pet has a low IgE level against a particular allergen. This can easily result in a false negative. In addition, extracts can vary not just between laboratories but also allergenic extract, making results hard to reproduce.

THE SOLUTION

To gain more accurate and sensitive information, tests are needed to identify each individual allergenic protein. Instead of testing for the house dust mite (or any particular allergen) as a whole, techniques are needed to test for the specific proteins that cause an allergic reaction. Human practitioners use molecular allergology to determine the allergens causing problems to provide their patients with a better level of care.

Macro Array Diagnostics launched the Allergy Explorer (ALEX), which provides a sensitization profile for human patients based on a test panel composed of allergen extracts and molecular allergens. Since their founding in 2016, the company has launched two generations of ALEX, offering a panel that covers nearly 100% of the world's relevant allergens. It also developed the Food Xplorer (FOX), to detect IgG-mediated food intolerances.

THE INNOVATION

After experiencing allergies myself, and being tested using molecular allergology, I thought the technology could be extremely beneficial to veterinary medicine.

This led to Nextmune partnering with Macro Array Diagnostics to develop the Pet Allergy Xplorer (PAX), the first commercial serological IgE-specific test that uses allergen extracts and molecular components to identify which allergens are affecting pets.

Advantages include:

- Improved reproducibility — A state-of-the-art robot builds each ELISA testing array, providing a the ability to produce an effective hyposensitization therapy that will provide more relief for the pet.
- Improved cross-reactivity identification — Using allergen extracts and molecular components in one test will help identify allergenic cross-reactivity.
- Elucidated poly sensitization — When multiple allergens cause a reaction, PAX will help identify the primary offenders.
- Individualized results — The results provided to the veterinarian will be specific to the region where they live, as well as the pet's species.
- Improved accuracy — PAX uses uniform production method that exceeds current ELISA plate building reproducibility. In addition, a standardization process is used to make the allergenic extracts, improving reproducibility.
- Increased data — Typical allergy testing via serum provides about 90 results, and intradermal testing provides around 60 to 80 results. The PAX cartridge holds 300 positions—100 of which will be extracts and 200 will be molecular components.
- Automated technology — Robots are used throughout the process, decreasing human error. The microarray dots during the plate manufacturing process are dispensed by a robot. In addition, a robot pipettes the serum during the testing phase. Controls are built into the PAX cartridge to ensure the sample has sufficient IgE for testing, and these controls also indicate if the cartridge is used properly.
- Improved treatment — The more accurate and sensitive results that PAX provides will improve pertinent allergens are identified to prevent confusion.
- Ongoing development — The cartridges currently contain the allergenic components and extracts that prevailing research indicates are clinically relevant. As data is gathered through testing, other allergenic components may be discovered. For example, if an extract continues to indicate positive while the individual molecular components indicate negative, we need to determine if another unidentified component is present or a cross reactivity is occurring. This data will help us characterize molecules to add in a single well-characterized anti-IgE monoclonal antibody to detect pet IgE, ensuring a detection level. PAX also uses technology to block cross-reactive carbohydrate determinants (CCD). In some cases, CCDs bind to IgE receptors, creating false positives. The PAX technology uses a blocking agent to prevent this binding. Allergy tests that do not use CCD blockers identify numerous allergens, many of which are irrelevant. This technology ensures only the future versions, revolutionizing the research in pet allergies.

Nextmune is the only veterinary diagnostic laboratory currently using molecular allergology. This next generation test will allow veterinarians to more accurately and sensitively diagnose allergic pets to facilitate treatment strategies. The initial launch for dogs is scheduled for January 2023, followed by tests for cats and horses.

Thierry Olivry, DrVet, PhD, DipECVD, DipACVD, Scientific Advisor and Project Leader
Dr. Olivry is a graduate of the University of Toulouse, France. He completed a dermatology residency and PhD in comparative pathology at UC Davis and is a Diplomate of both ACVD and ECVD. Dr. Olivry spent most of his career as a clinician-scientist at NC State University researching allergic and autoimmune skin diseases and now lives in Riga, Latvia.

 **nextmune**

